



**Knee care Plus**



**Knee Care Plus especially designed for knee with upgraded adhesive pads.**

■ Low Frequency Home Medical Device

Knee Care's unique supportive shape surrounds the knee joint and stimulates the area directly for improvement of knee pains and is the first home use treatment device specially designed for knees.


- Thermal therapy that softens tiredness, muscular pain, and neuralgia, and makes working of stomach and intestines active
- It is a sweep modulation wave current and a therapeutic gain to correspondence to a wide pain and a deeper part.
- Therapeutic instrument that wraps knee joint with supporter shape, and stimulates knee diseased part directly only for knee
- Moist sponge pad
- Will not miss any complex pains in knees

Product Name : Knee care Plus (Model name : SM2MT)

Japan Ministry of Health, Welfare, and Labor Certification Number : 222AGBZX00161A01


# New Sweeping & Modulating Waves of Electrical Current Technology

## Frequency




Therapeutic time


(Sweep) Gently changes the frequency of the waves.



(Modulation)  
Finely changes the size of the 4,000Hz waves and creates new waves of over frequencies (50 - 1000Hz).



## Electrode



The mechanisms for knee treatment are complicated involving various different frequencies.

The variable sweep currents give the therapeutic current at 4,000Hz AM variation to enhance the therapeutic effect and are auto-variable (sweep) between 50 and 1,000Hz.

This is our original technology which responds to a broad range of pain, and enriches the curative effect to a deeper part.

## Effects & Results

Loosens stiff shoulders, massages and prevents atrophy of numb muscles

(Massage Results:)

- Eases and soothes muscle pain and nerve pain
- Increases blood circulation
- Removes muscle fatigue
- Loosens stiff muscles

### ● Exclusive Knee Supporter

Electric therapy devices up until now has worked to improve blood flow in the thigh muscles in the lap, but not directly in the knees.

Knee Care Plus uses our original supporter that wraps around the knee joint, applies direct stimulation to the affect area, and relieves knee pain.



It's said that it is essential to do moderate exercise daily for easing pain and stiffness.

We highly recommend Knee Care Plus, an electric therapeutic device that is easy to use and provides comfortable stimulation.

- Place the electrode pads on your neck, back, and other places

Knee Care Plus can be used in various locations on the body. Effectively treat pain in places such as your neck, shoulders, back, elbows, and ankles through the electrode pads.